

# Biopor® AB Simply the best: When 2 become 1




Dear customers,

We would like to inform you that we will soon be able to offer you the best Biopor® AB ever. When 2 become 1: We have successfully combined the advantages of Biopor® AB Xtreme with those of Biopor® AB. Hence, this silicone offers you the best features of both products.




- Faster due to reduced curing time**
- Highest tensile strength**
- Wide color range**
- Curing at room temperature for easy handling**
- Highly precise, robust but also soft earmolds**
- Biocompatible**


Elongation at break, %



Tensile strength, MPa



Tear strength, N/mm



■ Biopor® AB (2024) – 40 Shore A   ■ Biopor® AB Xtreme – 40 Shore A   ■ Biopor® AB – 40 Shore A

Save valuable time and effort when processing the material: Curing time is reduced to 15 minutes and post-curing is possible at room temperature now! In addition, the colors are brighter and even more intense. The new Biopor® AB will still be available in the usual final hardnesses – 25, 40, 60 Shore A.

The conversion of the products is planned successively. Please do not hesitate to contact us if you have any questions or in case you are interested in sample material for test purposes – we will be happy to support you!

Kind regards

Your Dreve Otoplastik Team